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
The screenshot shows the TEM System Task Bar interface with the following components and callouts:


- HT/Beam Condition:** Acceleration voltage (160.00 kV) and Beam Current (81 uA). Callouts: 加速電圧On/OFF, ビームOn/OFF, 加速電圧.
- Illumination System:** Mode (TEM/SCAN), Spot Size (2), Alpha (2), and Illumination (TEM). Callouts: スポットサイズ, 収束角:  $\alpha$ , 照射モード (TEM/EDS/NBD/CBD), TEM/STEM切替.
- Stage:** Sample position (X, Y, Z, TX, TY) and speed controls. Callouts: 試料位置 (XY), 試料位置 (super fine X/Y), 試料高さ (Z), 試料傾斜 (TX/TY), スピード.
- Image Forming System:** Image mode (MAG1), Magnification (X50k), Defocus (0.0 nm). Callouts: 結像モード (MAG/DIFF), 倍率/カメラ長, デフォーカス量.
- Screen:** Current density (0.0 pA/cm<sup>2</sup>).

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


The screenshot displays the 'Condition Monitor' interface with three tabs: 'Status', 'Warning', and 'Log'. The 'Status' tab is active, showing a list of system components and their status. Each component has a black square checkbox to its left, which is currently unchecked. The components are grouped into four sections: 'Beam Control', 'Stage', 'Vacuum System', and 'Electron Optical System'. At the bottom right, there is a button labeled 'Alarm Stop'.

Section	Component	Status
Beam Control	Filament Over Current	Off
	High Voltage Over Current	Off
	High Voltage Cable Failure	Off
	Gun Gas Pressure Down	Off
Stage	X Motor Failure	Off
	Y Motor Failure	Off
	Z Motor Failure	Off
	Tilt X Motor Failure	Off
	Tilt Y / Rotation Motor Failure	Off
Vacuum System	Rotary Pump 1 Failure	Off
	Diffusion Pump 1 Disconnection	Off
	Diffusion Pump 1 Over Heat	Off
	Air Pressure Down	Off
	Reservoir Tank Vacuum Down	Off
	Pirani Gauge Disconnection	Off
Electron Optical System	OL Over Heat	Off
	OL Thermo Wire Break	Off

Alarm Stop

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- Beam
- Illumination System
- Specimen
- Stage
- Nanospace Map
- Aperture
- Image Forming System
- Screen
- Alignment Panel
- Wobbler

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

The screenshot shows the 'Beam Controller for JEM Administrators' interface. It features several control panels:

- HT Voltage:** Includes 'ON' and 'OFF' buttons, a 'Current HT' display at 160.00 [kV], and a row of preset buttons (80, 100, 120, 160, 200). The '160' button is highlighted.
- HT Voltage Settings:** A 'Target' field set to 160.00 [kV] and a 'Step' dropdown set to 1.0 [kV].
- Filament:** Includes 'ON' and 'OFF' buttons, and a 'Beam Current' display at 80.80 [uA].
- Filament Settings:** A 'Target' field set to 70.0 [%] and a 'Max' field set to 74 [%]. A 'Step' dropdown is set to 1.0 [%]. A progress bar shows the current filament current level.
- Filament Control Speed:** Includes 'Normal' and 'Degas' buttons, and a speed setting of 90 min.
- Bias:** Includes 'Bias Coarse' and 'Bias Fine' settings, both currently set to 7.

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
The image shows a software dialog box titled "Specimen". It contains a "Specimen Holder" section with a dropdown menu currently showing "EM-21010/21020 : Single Tilt Holder". Below this is an "Information" section with two text input fields: "Specimen:" and "Comment:". At the bottom of the dialog are two buttons: "OK" and "Cancel".

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

Alignment Panel for JEM Administrators

**Knob Assign**  
SHIFT-X,Y Knob : CL Alignment Shift  
DEF/STIG-X,Y Knob : (free)

**Deflector**  
GUN Align Spot Align Image Shift 1 Image Shift 2 PL Align  
Detect Align Dark Tilt Bright Tilt

**Stigmator**  
COND STIG OBJ STIG IL STIG

**Compensator**  
Shift Tilt Angle

**Neutral呼び出し**  
Neutralize Selected Alignment Neutralize All Alignment

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

The screenshot shows the 'HT Scheduling' control panel. It features a 'Linear Step' section with 'Step Volt' set to 0.1 [kV], 'Interval Time' set to 3 [s], and 'Total Time' set to 20.1 [min]. Below this, the 'Current HT Voltage' is 160.00 [kV] and the 'Target HT Voltage' is 200.00 [kV]. A 'Start' button is highlighted with a red box and labeled '昇圧開始' (Start Voltage Increase). A 'Stop' button is also present. A red box with a warning symbol and 'HT OFF' is labeled '加速電圧OFF ->0kV' (Acceleration Voltage OFF -> 0kV). At the bottom, an 'HT Conditioning' section has an 'ON' button and a 'Hold' checkbox, both highlighted with a red box and labeled '管理者用' (Administrator Use).

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